


## Kaydee Sterile Disposable Skin Stapler and Staple Remover



\*Angled head provides clear view to ensure precise staple placement for staple to enter tissue easily

\*Quick staple release

\*Easy to use/lightweight /handle designed for small to large hands

\* Suggested uses - General Surgery for wound closure

### Sterile Disposable Skin Stapler



Two skin stapler models are Model D and E, Three staple sizes

- 1) *Wide* - general part of body - 7cm wide x 4cm high x 0.6cm staple mouth
- 2) *Regular* - thin part of body ie Face/Head - 5.4cm wide x 3.6 high x 0.5cm staple mouth
- 3) *High* -parts of body such as Stomach - 7cm wide x 4.5cm high x 0.6cm staple mouth

- Material of staple – 316L stainless steel or titanium (better for MRI)
- Nine stapler counts 15/25/35/20/30/40/45/50/55
- Sterile to 3 years from date of Manufacture
- For single use only
- Packaging 80 pieces for carton
- Countries that use the Kaydee Stapler and Staple remover include USA, Turkey, Korea, India, Hong Kong, Italy, France, Australia
- Designed in Germany and manufactured in China
- TGA Approved – ARTG Number – 209117

## Kaydee Sterile Disposable Staple Remover



- Can remove all surgical skin staple types- all models do the same role
- Three model Staple Remover
  - a) T-1 Model – single use
  - b) T-2 Model - single use – most commonly used
  - c) T-3 Model stainless steel and reusable
- Sterile to 3 years from date of Manufacture
- Packaging 300 pieces per carton

*Manufactured by Xian Kaydee Medical Appliances Co Ltd  
8/F Building 3, No 1 Torch Road, Xian 710043, Shaanxi, China*

*Distributed by Global Surgical Distribution  
115 The Promenade Camp Hill, Queensland, Australia 4152  
Phone: 61 7 39016311  
Mobile: 040 1969 615  
ABN: 15 152 255 976*

